

The "B" Clip (Bifurcation Clip)

Published date: Oct. 23, 2013

Technology description

One of the complications of EVAR devices used to treat abdominal aortic aneurysms endovascularly is distal migration. Approximation of the iliac limbs of the device using a special clip applied during the EVAR implant procedure, or as a subsequent procedure, will help to prevent or halt distal device migration, which is a cause of device failure.

Institution

[University of Oklahoma](#)

Inventors

[Jorge Saucedo](#)

Professor

Medicine/Cardiology

[Wael Abo-Auda](#)

Resident

Dept. of Medicine Cardiology

联系我们



叶先生

电话 : 021-65679356

手机 : 13414935137

邮箱 : yeyingsheng@zf-ym.com